

Mcl-1 (phospho-Ser64) rabbit pAb

Catalog No: #13728



Package Size: #13728-1 50ul #13728-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Mcl-1 (phospho-Ser64) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human Mcl-1 (phospho-Ser64)
Immunogen Description	Synthesized phosho peptide around human Mcl-1 (Ser64)
Other Names	Induced myeloid leukemia cell differentiation protein Mcl-1 (Bcl-2-like protein 3) (Bcl2-L-3) (Bcl-2-related protein EAT/mcl1) (mcl1/EAT)
Accession No.	Swiss Prot:Q07820GeneID:4170
Uniprot	Q07820
GeneID	4170
SDS-PAGE MW	37
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Background

BCL2 family apoptosis regulator(MCL1) Homo sapiens This gene encodes an anti-apoptotic protein, which is a member of the Bcl-2 family. Alternative splicing results in multiple transcript variants. The longest gene product (isoform 1) enhances cell survival by inhibiting apoptosis while the alternatively spliced shorter gene products (isoform 2 and isoform 3) promote apoptosis and are death-inducing. [provided by RefSeq, Oct 2010],

Note: This product is for in vitro research use only