

Lyn (phospho-Tyr507) rabbit pAb

Catalog No: #13740

Package Size: #13740-1 50ul #13740-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lyn (phospho-Tyr507) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human,Mouse
Specificity	This antibody detects endogenous levels of Human Mouse Lyn (phospho-Tyr507)
Immunogen Description	Synthesized phosho peptide around human Lyn (Tyr507)
Other Names	Tyrosine-protein kinase Lyn (EC 2.7.10.2) (Lck/Yes-related novel protein tyrosine kinase) (V-yes-1 Yamaguchi sarcoma viral related oncogene homolog) (p53Lyn) (p56Lyn)
Accession No.	Swiss Prot:P07948GenelD:4067
Uniprot	P07948
GeneID	4067
SDS-PAGE MW	58
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Background

LYN proto-oncogene, Src family tyrosine kinase(LYN) Homo sapiens This gene encodes a tyrosine protein kinase, which maybe involved in the regulation of mast cell degranulation, and erythroid differentiation. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2011],

Note: This product is for in vitro research use only