

IRF-7 (phospho-Ser471/472) rabbit pAb

Catalog No: #13785



Package Size: #13785-1 50ul #13785-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IRF-7 (phospho-Ser471/472) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human IRF-7 (phospho-Ser471 or 472)
Immunogen Description	Synthesized phospho peptide around human IRF-7 (Ser471 and 472)
Other Names	Interferon regulatory factor 7 (IRF-7)
Accession No.	Swiss Prot:Q92985GeneID:3665
Uniprot	Q92985
GeneID	3665
SDS-PAGE MW	55
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Background

interferon regulatory factor 7(IRF7) Homo sapiens IRF7 encodes interferon regulatory factor 7, a member of the interferon regulatory transcription factor (IRF) family. IRF7 has been shown to play a role in the transcriptional activation of virus-inducible cellular genes, including interferon beta chain genes. Inducible expression of IRF7 is largely restricted to lymphoid tissue. Multiple IRF7 transcript variants have been identified, although the functional consequences of these have not yet been established. [provided by RefSeq, Jul 2008],

Note: This product is for in vitro research use only