

GARNL1 (Phospho-S797) Polyclonal Antibody

Catalog No: #13843



Package Size: #13843-1 50ul #13843-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GARNL1 (Phospho-S797) Polyclonal Antibody
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	IHC-p,IF(paraffin section),WB
Species Reactivity	Human, Mouse, Rat
Specificity	This antibody detects endogenous phospho levels of GARNL1 (Phospho-S797) at Human, Mouse, Rat
Immunogen Description	Synthesized peptide derived from human GARNL1 (Phospho-S797)
Other Names	Ral GTPase-activating protein subunit alpha-1 (GAP-related-interacting partner to E12) (GRIPE) (GTPase-activating Rap/Ran-GAP domain-like 1) (Tuberin-like protein 1) (p240)
Accession No.	Swiss Prot:Q6GYQ0GenelD:253959
Uniprot	Q6GYQ0
GeneID	253959
SDS-PAGE MW	220
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

IHC-p 1:50-200, WB 1:500-2000

Background

Ral GTPase activating protein catalytic alpha subunit 1(RALGAPA1) Homo sapiens This gene encodes a major subunit of the RAL-GTPase activating protein. A similar protein in mouse binds E12, a transcriptional regulator of immunoglobulin genes. The mouse protein also functions in skeletal muscle by binding to the regulatory 14-3-3 proteins upon stimulation with insulin or muscle contraction. A pseudogene of this gene has been identified on chromosome 9. [provided by RefSeq, Oct 2016].

Note: This product is for in vitro research use only