

FRS2- α (phospho-Tyr196) rabbit pAb

Catalog No: #13845



Package Size: #13845-1 50ul #13845-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FRS2- α (phospho-Tyr196) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human, Mouse
Specificity	This antibody detects endogenous levels of Human Mouse FRS2- α (phospho-Tyr196)
Immunogen Description	Synthesized phospho peptide around human FRS2- α (Tyr196)
Other Names	Fibroblast growth factor receptor substrate 2 (FGFR substrate 2) (FGFR-signaling adaptor SNT) (Suc1-associated neurotrophic factor target 1) (SNT-1)
Accession No.	Swiss Prot:Q8WU20 GenID:10818
Uniprot	Q8WU20
GeneID	10818
SDS-PAGE MW	65
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only