

FRA1 (phospho-Ser265) rabbit pAb

Catalog No: #13846



Package Size: #13846-1 50ul #13846-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FRA1 (phospho-Ser265) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	This antibody detects endogenous levels of Human Mouse Rat FRA1 (phospho-Ser265)
Immunogen Description	Synthesized phospho peptide around human FRA1 (Ser265)
Other Names	Fos-related antigen 1 (FRA-1)
Accession No.	Swiss Prot:P15407GenelD:8061
Uniprot	P15407
GenelD	8061
SDS-PAGE MW	43
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Background

FOS like 1, AP-1 transcription factor subunit(FOSL1) Homo sapiens The Fos gene family consists of 4 members: FOS, FOSB, FOSL1, and FOSL2. These genes encode leucine zipper proteins that can dimerize with proteins of the JUN family, thereby forming the transcription factor complex AP-1. As such, the FOS proteins have been implicated as regulators of cell proliferation, differentiation, and transformation. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2014],

Note: This product is for in vitro research use only