

FoxL2 (phospho Ser263) Polyclonal Antibody

Catalog No: #13856



Package Size: #13856-1 50ul #13856-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FoxL2 (phospho Ser263) Polyclonal Antibody
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Specificity	Phospho-FoxL2 (S263) Polyclonal Antibody detects endogenous levels of FoxL2 protein only when phosphorylated at S263.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human FOXL2 around the phosphorylation site of Ser263. AA range:229-278
Other Names	FOXL2; Forkhead box protein L2
Accession No.	Swiss Prot:P58012GeneID:668
Uniprot	P58012
GeneID	668
SDS-PAGE MW	40
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.

Background

forkhead box L2(FOXL2) Homo sapiens This gene encodes a forkhead transcription factor. The protein contains a fork-head DNA-binding domain and may play a role in ovarian development and function. Expansion of a polyalanine repeat region and other mutations in this gene are a cause of blepharophimosis syndrome and premature ovarian failure 3. [provided by RefSeq, Jul 2016].

Note: This product is for in vitro research use only