

eIF4B (phospho-Ser406) rabbit pAb

Catalog No: #13900



Package Size: #13900-1 50ul #13900-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	eIF4B (phospho-Ser406) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	This antibody detects endogenous levels of Human Mouse Rat eIF4B (phospho-Ser406)
Immunogen Description	Synthesized phospho peptide around human eIF4B (Ser406)
Other Names	Eukaryotic translation initiation factor 4B (eIF-4B)
Accession No.	Swiss Prot:P23588GenelD:1975
Uniprot	P23588
GenelD	1975
SDS-PAGE MW	80
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only