

DR6 (phospho-Ser562) rabbit pAb

Catalog No: #13910



Package Size: #13910-1 50ul #13910-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DR6 (phospho-Ser562) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human,Mouse,Rat
Specificity	This antibody detects endogenous levels of Human Mouse Rat DR6 (phospho-Ser562)
Immunogen Description	Synthesized phosho peptide around human DR6 (Ser562)
Other Names	Tumor necrosis factor receptor superfamily member 21 (Death receptor 6) (CD antigen CD358)
Accession No.	Swiss Prot:O75509GeneID:27242
Uniprot	O75509
GeneID	27242
SDS-PAGE MW	72
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Background

TNF receptor superfamily member 21(TNFRSF21) Homo sapiens This gene encodes a member of the tumor necrosis factor receptor superfamily. The encoded protein activates nuclear factor kappa-B and mitogen-activated protein kinase 8 (also called c-Jun N-terminal kinase 1), and induces cell apoptosis. Through its death domain, the encoded receptor interacts with tumor necrosis factor receptor type 1-associated death domain (TRADD) protein, which is known to mediate signal transduction of tumor necrosis factor receptors. Knockout studies in mice suggest that this gene plays a role in T-helper cell activation, and may be involved in inflammation and immune regulation. [provided by RefSeq, Jul 2013],

Note: This product is for in vitro research use only