

DBC1 (phospho-Thr454) rabbit pAb

Catalog No: #13925

Package Size: #13925-1 50ul #13925-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DBC1 (phospho-Thr454) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human DBC1 (phospho-Thr454)
Immunogen Description	Synthesized phospho peptide around human DBC1 (Thr454)
Other Names	DBIRD complex subunit KIAA1967 (Deleted in breast cancer gene 1 protein) (DBC-1) (DBC.1) (p30 DBC)
Accession No.	Swiss Prot:Q8N163GeneID:57805
Uniprot	Q8N163
GeneID	57805
SDS-PAGE MW	105
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only