

DAPP1/BAM32 (phospho-Tyr139) rabbit pAb

Catalog No: #13927



Package Size: #13927-1 50ul #13927-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	DAPP1/BAM32 (phospho-Tyr139) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human DAPP1 or BAM32 (phospho-Tyr139)
Immunogen Description	Synthesized phosho peptide around human DAPP1 (Tyr139)
Other Names	Dual adapter for phosphotyrosine and 3-phosphotyrosine and 3-phosphoinositide (hDAPP1) (B lymphocyte adapter protein Bam32) (B-cell adapter molecule of 32 kDa)
Accession No.	Swiss Prot:Q9UN19GeneID:27071
Uniprot	Q9UN19
GeneID	27071
SDS-PAGE MW	32
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only