

CTDSPL2 (phospho-Ser104) rabbit pAb

Catalog No: #13938



Package Size: #13938-1 50ul #13938-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CTDSPL2 (phospho-Ser104) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human CTDSPL2 (phospho-Ser104)
Immunogen Description	Synthesized phospho peptide around human CTDSPL2 (Ser104)
Other Names	CTD small phosphatase-like protein 2 (CTDSP-like 2) (EC 3.1.3.-)
Accession No.	Swiss Prot:Q05D32GeneID:51496
Uniprot	Q05D32
GeneID	51496
SDS-PAGE MW	50
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only