

BCL9L (phospho-Ser915) rabbit pAb

Catalog No: #14050



Package Size: #14050-1 50ul #14050-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|-----------------------------------------------------------------------------------------------------------|
| Product Name | BCL9L (phospho-Ser915) rabbit pAb |
| Host Species | Rabbit |
| Purification | The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen. |
| Applications | WB |
| Species Reactivity | Human |
| Specificity | This antibody detects endogenous levels of Human BCL9L (phospho-Ser915) |
| Immunogen Description | Synthesized phospho peptide around human BCL9L (Ser915) |
| Other Names | B-cell CLL/lymphoma 9-like protein (B-cell lymphoma 9-like protein) (BCL9-like protein) (Protein BCL9-2) |
| Accession No. | Swiss Prot:Q86UU0GeneID:283149 |
| Uniprot | Q86UU0 |
| GeneID | 283149 |
| SDS-PAGE MW | 160 |
| Concentration | 1 mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | -20°C/1 |

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only