

## BACH1/BRIP1 (phospho-Thr1133) rabbit pAb

Catalog No: #14056



Package Size: #14056-1 50ul #14056-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	BACH1/BRIP1 (phospho-Thr1133) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human BACH1 or BRIP1 (phospho-Thr1133)
Immunogen Description	Synthesized phospho peptide around human BACH1 and BRIP1 (Thr1133)
Other Names	Fanconi anemia group J protein (Protein FACJ) (EC 3.6.4.13) (ATP-dependent RNA helicase BRIP1) (BRCA1-associated C-terminal helicase 1) (BRCA1-interacting protein C-terminal helicase 1) (BRCA1-interacting protein 1)
Accession No.	Swiss Prot:Q9BX63GeneID:83990
Uniprot	Q9BX63
GeneID	83990
SDS-PAGE MW	135
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

## Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only