

Atg14 (phospho-Ser29) rabbit pAb

Catalog No: #14061



Package Size: #14061-1 50ul #14061-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Atg14 (phospho-Ser29) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human,Mouse
Specificity	This antibody detects endogenous levels of Human Mouse Atg14 (phospho-Ser29)
Immunogen Description	Synthesized phospho peptide around human Atg14 (Ser29)
Other Names	Beclin 1-associated autophagy-related key regulator (Barkor) (Autophagy-related protein 14-like protein) (Atg14L)
Accession No.	Swiss Prot:Q6ZNE5GeneID:22863
Uniprot	Q6ZNE5
GeneID	22863
SDS-PAGE MW	55
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only