

## Atg13 (phospho-Ser355) rabbit pAb

Catalog No: #14062

Package Size: #14062-1 50ul #14062-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Atg13 (phospho-Ser355) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of Human Atg13 (phospho-Ser355)
Immunogen Description	Synthesized phospho peptide around human Atg13 (Ser355)
Other Names	Autophagy-related protein 13
Accession No.	Swiss Prot:O75143GeneID:9776
Uniprot	O75143
GeneID	9776
SDS-PAGE MW	57
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

## Application Details

WB 1:1000-2000

## Background

autophagy related 13(ATG13) Homo sapiens The protein encoded by this gene is an autophagy factor and a target of the TOR kinase signaling pathway. The encoded protein is essential for autophagosome formation and mitophagy. [provided by RefSeq, Oct 2016],

Note: This product is for in vitro research use only