

53BP1 (phospho-Ser1618) rabbit pAb

Catalog No: #14115

Package Size: #14115-1 50ul #14115-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	53BP1 (phospho-Ser1618) rabbit pAb
Host Species	Rabbit
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB
Species Reactivity	Human,Rat
Specificity	This antibody detects endogenous levels of Human Rat 53BP1 (phospho-Ser1618)
Immunogen Description	Synthesized phospho peptide around human 53BP1 (Ser1618)
Other Names	Tumor suppressor p53-binding protein 1 (53BP1) (p53-binding protein 1) (p53BP1)
Accession No.	Swiss Prot:Q12888GeneID:7158
Uniprot	Q12888
GeneID	7158
SDS-PAGE MW	213
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	-20°C/1

Application Details

WB 1:1000-2000

Note: This product is for in vitro research use only