

# Recombinant Human SH3 domain and tetratricopeptide repeat-containing protein 1(SH3TC1),partial

Catalog No: #AP0102

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP0102-1 20ug #AP0102-2 100ug #AP0102-3 1mg

## Description

Product Name	Recombinant Human SH3 domain and tetratricopeptide repeat-containing protein 1(SH3TC1),partial
Brief Description	Recombinant Proteins
Host Species	Yeast
Target Name	His-tagged
Other Names	SH3TC1SH3 domain and tetratricopeptide repeat-containing protein 1
Tag Info	N-terminal His-tagged
Target Sequence	1-482aa
Concentration	0.1-1.0 mg/mL,by the Bradford Method.
Formulation	Tris-based buffer,50% glycerol
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>

## Application Details

Greater than 85% as determined by SDS-PAGE.

Note: This product is for in vitro research use only