

Recombinant Human Fibronectin type III domain-containing protein 3B(FNDC3B) ,partial

Catalog No: #AP0103

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP0103-1 20ug #AP0103-2 100ug #AP0103-3 1mg

Description

| | |
|-------------------|--|
| Product Name | Recombinant Human Fibronectin type III domain-containing protein 3B(FNDC3B) ,partial |
| Brief Description | Recombinant Proteins |
| Host Species | Yeast |
| Target Name | His-tagged |
| Other Names | FNDC3B; FAD104; NS5ABP37; UNQ2421/PRO4979/PRO34274Fibronectin type III domain-containing protein 3B; Factor for adipocyte differentiation 104; HCV NS5A-binding protein 37 |
| Tag Info | N-terminal His-tagged |
| Target Sequence | 278-570aa |
| Concentration | 0.1-1.0 mg/mL,by the Bradford Method. |
| Formulation | Tris-based buffer,50% glycerol |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |

Application Details

Greater than 85% as determined by SDS-PAGE.

Background

May be a positive regulator of adipogenesis.

Note: This product is for in vitro research use only