

Recombinant Human Transcription factor GATA-6(GATA6)

Catalog No: #AP0106

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP0106-1 20ug #AP0106-2 100ug #AP0106-3 1mg

Description

Product Name	Recombinant Human Transcription factor GATA-6(GATA6)
Brief Description	Recombinant Proteins
Host Species	Yeast
Target Name	His-tagged
Other Names	Gata binding factor 6; Gata binding protein 6; GATA-binding factor 6; Gata6; GATA6_HUMAN; Transcription factor Gata 6; Transcription factor GATA-6
Tag Info	N-terminal His-tagged
Target Sequence	1-595aa
Concentration	0.1-1.0 mg/mL,by the Bradford Method.
Formulation	Tris-based buffer,50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Application Details

Greater than 85% as determined by SDS-PAGE.

Background

Transcriptional activator

Note: This product is for in vitro research use only