

# Recombinant Human Transcription factor SOX-9(SOX9)

Catalog No: #AP0108

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP0108-1 20ug #AP0108-2 100ug #AP0108-3 1mg

## Description

Product Name	Recombinant Human Transcription factor SOX-9(SOX9)
Brief Description	Recombinant Proteins
Host Species	Yeast
Target Name	His-tagged
Other Names	campomelic dysplasia autosomal sex reversal; CMD 1; CMD1; CMPD 1; CMPD1; SOX 9; Sox9; SOX9_HUMAN; SRA 1; SRA1; SRXX2; SRXY10; SRY (sex determining region Y) box 9 (campomelic dysplasia autosomal; SRY (sex determining region Y) box 9; SRY (sex determining region Y)-box 9; SRY (sex-determining region Y)-box 9 protein; SRY related HMG box gene 9; Transcription factor SOX 9; Transcription factor SOX-9; transcription factor SOX9
Tag Info	N-terminal His-tagged
Target Sequence	1-509aa
Concentration	0.1-1.0 mg/mL,by the Bradford Method.
Formulation	Tris-based buffer,50% glycerol
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.  Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

## Application Details

Greater than 85% as determined by SDS-PAGE.

## Background

Transcriptional regulator

Note: This product is for in vitro research use only