

# Human GNPAT1 / Glucosamine 6-phosphate N-acetyltransferase ELISA Kit



Catalog No: #EK21001

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human GNPAT1 / Glucosamine 6-phosphate N-acetyltransferase ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Phosphoglucosamine acetylase, Phosphoglucosamine transacetylase, GNA1
Accession No.	Uniprot:Q96EK6
Uniprot	Q96EK6
Storage	6 months

## Application Details

Detection Range: 0.312-20ng/mL Assay Type: Sandwich Sensitivity: 0.101ng/mL

## Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only