

# Human DAAM1 / Disheveled-associated activator of morphogenesis 1 ELISA Kit



Catalog No: #EK21002

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human DAAM1 / Disheveled-associated activator of morphogenesis 1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	KIAA0666
Accession No.	Uniprot: Q9Y4D1
Uniprot	Q9Y4D1
Storage	6 months

## Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.098ng/mL

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only