

Human FAM83B / Protein FAM83B ELISA Kit

Catalog No: #EK21003



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human FAM83B / Protein FAM83B ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	C6orf143
Accession No.	Uniprot:Q5T0W9
Uniprot	Q5T0W9
Storage	6 months

Application Details

Detection Range:0.78-50ng/mL Assay Type:Sandwich Sensitivity:0.24ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only