

# Human DUOX1 / Dual oxidase 1 ELISA Kit

Catalog No: #EK21004



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Human DUOX1 / Dual oxidase 1 ELISA Kit  |
| Brief Description  | ELISA Kit   |
| Applications       | ELISA   |
| Species Reactivity | Human   |
| Other Names        | Large NOX 1, Long NOX 1, NADPH thyroid oxidase 1, Thyroid oxidase 1, DUOX, LNOX1, THOX1 |
| Accession No.      | Uniprot:Q9NRD9  |
| Uniprot            | Q9NRD9  |
| Storage            | 6 months  |

## Application Details

Detection Range: 0.156-10ng/mL Assay Type: Sandwich Sensitivity: 0.066ng/mL

## Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only