

Human VASH2 / Tubuliny-Tyr carboxypeptidase 2 ELISA Kit



Catalog No: #EK21006

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human VASH2 / Tubuliny-Tyr carboxypeptidase 2 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	VASHL, Vasohibin-2, Vasohibin-like protein
Accession No.	Uniprot:Q86V25
Uniprot	Q86V25
Storage	6 months

Application Details

Detection Range:0.78-50ng/mL Assay Type:Sandwich Sensitivity:0.24ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only