

# Rat Hfe / Hereditary hemochromatosis protein homolog ELISA Kit



Catalog No: #EK21012

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Rat Hfe / Hereditary hemochromatosis protein homolog ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	RT
Other Names	RT1-CAFE
Accession No.	Uniprot:O35799
Uniprot	O35799
Storage	6 months

## Application Details

Detection Range:0.78-50ng/mL Assay Type:Sandwich Sensitivity:0.27ng/mL

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only