

Human ILF2 / Interleukin enhancer-binding factor 2 ELISA Kit



Catalog No: #EK21013

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human ILF2 / Interleukin enhancer-binding factor 2 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Nuclear factor of activated T-cells 45 kDa,NF45,PRO3063
Accession No.	Uniprot:Q12905
Uniprot	Q12905
Storage	6 months

Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.085ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only