

Human ACMSD / 2-amino-3-carboxymuconate-6-semialdehyde decarboxylase ELISA Kit

Catalog No: #EK21017

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human ACMSD / 2-amino-3-carboxymuconate-6-semialdehyde decarboxylase ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Picolinate carboxylase
Accession No.	Uniprot:Q8TDX5
Uniprot	Q8TDX5
Storage	6 months

Application Details

Detection Range:0.78-50ng/mL Assay Type:Sandwich Sensitivity:0.29ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only