

# Human CNTNAP4 / Contactin-associated protein-like 4 ELISA Kit



Catalog No: #EK21018

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human CNTNAP4 / Contactin-associated protein-like 4 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Cell recognition molecule Caspr4,CASPR4,KIAA1763
Accession No.	Uniprot:Q9C0A0
Uniprot	Q9C0A0
Storage	6 months

## Application Details

Detection Range:0.156-10ng/mL Assay Type:Sandwich Sensitivity:0.056ng/mL

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only