

Mouse Mtap / S-methyl-5'-thioadenosine phosphorylase ELISA Kit



Catalog No: #EK21019

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Mtap / S-methyl-5'-thioadenosine phosphorylase ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse
Other Names	5'-methylthioadenosine phosphorylase,MTA phosphorylase
Accession No.	Uniprot:Q9CQ65
Uniprot	Q9CQ65
Storage	6 months

Application Details

Detection Range:78-5000pg/mL Assay Type:Sandwich Sensitivity:39pg/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only