

Human CSRP3 / Cysteine and glycine-rich protein 3 ELISA Kit



Catalog No: #EK21020

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human CSRP3 / Cysteine and glycine-rich protein 3 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Cardiac LIM protein, Cysteine-rich protein 3, CRP3, LIM domain protein, cardiac, Muscle LIM protein, CLP, MLP
Accession No.	Uniprot:P50461
Uniprot	P50461
Storage	6 months

Application Details

Detection Range: 0.156-10ng/mL Assay Type: Sandwich Sensitivity: 0.078ng/mL

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only