

Human PEAK1 / Inactive tyrosine-protein kinase PEAK1 ELISA Kit



Catalog No: #EK21024

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human PEAK1 / Inactive tyrosine-protein kinase PEAK1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	KIAA2002,Pseudopodium-enriched atypical kinase 1,Sugen kinase 269,Tyrosine-protein kinase SgK269
Accession No.	Uniprot:Q9H792
Uniprot	Q9H792
Storage	6 months

Application Details

Detection Range:78-5000pg/mL Assay Type:Sandwich Sensitivity:24pg/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only