

Human GDF7 / Growth/differentiation factor 7 ELISA Kit



Catalog No: #EK21025

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human GDF7 / Growth/differentiation factor 7 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	GDF-7
Accession No.	Uniprot:Q7Z4P5
Uniprot	Q7Z4P5
Storage	6 months

Application Details

Detection Range:78-5000pg/mL Assay Type:Sandwich Sensitivity:23pg/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only