

Human CFHR5 / Complement factor H-related protein 5 ELISA Kit



Catalog No: #EK21026

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human CFHR5 / Complement factor H-related protein 5 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	FHR-5,CFHL5,FHR5
Accession No.	Uniprot:Q9BXR6
Uniprot	Q9BXR6
Storage	6 months

Application Details

Detection Range:78-5000pg/mL Assay Type:Sandwich Sensitivity:30pg/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only