

# Human LFNG / Beta-1,3-N-acetylglucosaminyltransferase lunatic fringe ELISA Kit

Catalog No: #EK21029

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human LFNG / Beta-1,3-N-acetylglucosaminyltransferase lunatic fringe ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	O-fucosylpeptide 3-beta-N-acetylglucosaminyltransferase
Accession No.	Uniprot:Q8NES3
Uniprot	Q8NES3
Storage	6 months

## Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.15ng/mL

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only