

Mouse Nlgn1 / Neuroligin-1 ELISA Kit

Catalog No: #EK21033



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Nlgn1 / Neuroligin-1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse
Other Names	Kiaa1070
Accession No.	Uniprot:Q99K10
Uniprot	Q99K10
Storage	6 months

Application Details

Detection Range:0.156-10ng/mL Assay Type:Sandwich Sensitivity:0.078ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only