

Human UBXN7 / UBX domain-containing protein 7 ELISA Kit



Catalog No: #EK21034

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human UBXN7 / UBX domain-containing protein 7 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	KIAA0794,UBXD7
Accession No.	Uniprot:O94888
Uniprot	O94888
Storage	6 months

Application Details

Detection Range:0.78-50ng/mL Assay Type:Sandwich Sensitivity:0.41ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only