

# Human GKN1 / Gastroke-1 ELISA Kit

Catalog No: #EK21035



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human GKN1 / Gastroke-1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	18 kDa antrum mucosa protein,AMP-18,Protein CA11,AMP18,CA11,UNQ489/PRO1005
Accession No.	Uniprot:Q9NS71
Uniprot	Q9NS71
Storage	6 months

## Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.186ng/mL

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only