

Human CCL1 / C-C motif chemokine 1 ELISA Kit

Catalog No: #EK21037



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human CCL1 / C-C motif chemokine 1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Small-inducible cytokine A1,T lymphocyte-secreted protein I-309,SCYA1
Accession No.	Uniprot:P22362
Uniprot	P22362
Storage	6 months

Application Details

Detection Range:78-5000pg/mL Assay Type:Sandwich Sensitivity:39pg/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only