

Human NPC1 / NPC intracellular cholesterol transporter 1 ELISA Kit



Catalog No: #EK21038

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human NPC1 / NPC intracellular cholesterol transporter 1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Niemann-Pick C1 protein
Accession No.	Uniprot:O15118
Uniprot	O15118
Storage	6 months

Application Details

Detection Range:78-5000pg/mL Assay Type:Sandwich Sensitivity:39pg/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only