

Human ONECUT1 / Hepatocyte nuclear factor 6 ELISA Kit



Catalog No: #EK21039

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human ONECUT1 / Hepatocyte nuclear factor 6 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	HNF-6,One cut domain family member 1,One cut homeobox 1,HNF6,HNF6A
Accession No.	Uniprot:Q9UBC0
Uniprot	Q9UBC0
Storage	6 months

Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.078ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only