

Human DIAPH1 / Protein diaphanous homolog 1 ELISA Kit



Catalog No: #EK21041

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human DIAPH1 / Protein diaphanous homolog 1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Diaphanous-related formin-1,DRF1,DIAP1
Accession No.	Uniprot:O60610
Uniprot	O60610
Storage	6 months

Application Details

Detection Range:0.312-20ng/mL Assay Type:Sandwich Sensitivity:0.098ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only