

Human CLOCK / Circadian locomoter output cycles protein kaput ELISA Kit



Catalog No: #EK21042

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human CLOCK / Circadian locomoter output cycles protein kaput ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	hCLOCK,Class E basic helix-loop-helix protein 8,bHLHe8,BHLHE8,KIAA0334
Accession No.	Uniprot:O15516
Uniprot	O15516
Storage	6 months

Application Details

Detection Range:0.156-10ng/mL Assay Type:Sandwich Sensitivity:0.089ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only