

# Human AQP7 / Aquaporin-7 ELISA Kit

Catalog No: #EK21043



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human AQP7 / Aquaporin-7 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	AQP-7, Aquaglyceroporin-7, Aquaporin adipose, AQPap, Aquaporin-7-like, AQP7L, AQP9
Accession No.	Uniprot:O14520
Uniprot	O14520
Storage	6 months

## Application Details

Detection Range: 0.312-20ng/mL Assay Type: Sandwich Sensitivity: 0.183ng/mL

## Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only