

Human VAT1 / Synaptic vesicle membrane protein VAT-1 homolog ELISA Kit



Catalog No: #EK21044

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human VAT1 / Synaptic vesicle membrane protein VAT-1 homolog ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Accession No.	Uniprot:Q99536
Uniprot	Q99536
Storage	6 months

Application Details

Detection Range:0.156-10ng/mL Assay Type:Sandwich Sensitivity:0.078ng/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only