

Human FAM72A / Protein FAM72A ELISA Kit

Catalog No: #EK21045



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human FAM72A / Protein FAM72A ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Latent membrane protein 1-induced protein,LMP1-induced protein,UGENE
Accession No.	Uniprot:Q5TYM5
Uniprot	Q5TYM5
Storage	6 months

Application Details

Detection Range:78-5000pg/mL Assay Type:Sandwich Sensitivity:38pg/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only