

Mouse Chrac1 / Chromatin accessibility complex protein 1 ELISA Kit



Catalog No: #EK21047

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Chrac1 / Chromatin accessibility complex protein 1 ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse
Other Names	CHRAC-1, DNA polymerase epsilon subunit p15, NF-YC-like protein, YC-like protein 1, YCL1
Accession No.	Uniprot:Q9JKP8
Uniprot	Q9JKP8
Storage	6 months

Application Details

Detection Range: 7.8-500pg/mL Assay Type: Sandwich Sensitivity: 3.82pg/mL

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only