

Human HLF / Hepatic leukemia factor ELISA Kit

Catalog No: #EK21048



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human HLF / Hepatic leukemia factor ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Accession No.	Uniprot:Q16534
Uniprot	Q16534
Storage	6 months

Application Details

Detection Range:78-5000pg/mL Assay Type:Sandwich Sensitivity:35pg/mL

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only