

Human OPTN / Optineurin ELISA Kit

Catalog No: #EK21049

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human OPTN / Optineurin ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	E3-14.7K-interacting protein, FIP-2, Huntingtin yeast partner L, Huntingtin-interacting protein 7, HIP-7, Huntingtin-interacting protein L, NEMO-related protein, Optic neuropathy-inducing protein, Transcription factor IIIA-interacting protein, TFIIIA-IntP, FIP2, GLC1E, HIP7, HYPL, NRP
Accession No.	Uniprot:Q96CV9
Uniprot	Q96CV9
Storage	6 months

Application Details

Detection Range: 0.312-20ng/mL Assay Type: Sandwich Sensitivity: 0.163ng/mL

Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only